

## Biotinylated Anti-cynomolgus monkey CXCL9 Polyclonal Antibody

**KPAB0358**

### Description

---

This Biotinylated Anti-cynomolgus monkey CXCL9 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

### Product Information

---

**SKU:** KPAB0358

**Conjugated:** Biotin

**Applications:** ELISA, ELISpot

### Antibody Data

---

**Host:** Rabbit

**Reactivity:** Cynomolgus Monkey

**Immunogen:** Recombinant Cynomolgus Monkey CXCL9

**Purification:** Antigen-affinity Purification

**Homology:** Chlorocebus sabaeus (green monkey), Macaca fascicularis (crab-eating macaque), Macaca nemestrina (pig-tailed macaque), Mandrillus leucophaeus (drill), Theropithecus gelada (gelada)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.