

Anti-FNDC5 Polyclonal Antibody

KPAB0363

Description

This Anti-FNDC5 Polyclonal Antibody has been validated for WB, ELISA applications.

Product Information

SKU: KPAB0363
Conjugated: Unconjugated
Applications: WB, ELISA

Antibody Data

Host: Rabbit
Reactivity: Multi-species
Immunogen: Recombinant FNDC5
Purification: Antigen-affinity Purification

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.