

Biotinylated Anti-FNDC5 Polyclonal Antibody

KPAB0364

Description

This Biotinylated Anti-FNDC5 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0364
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit
Reactivity: Multi-species
Immunogen: Recombinant FNDC5
Purification: Antigen-affinity Purification

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.