

Anti-mouse IFN Alpha Polyclonal Antibody

KPAB0369

Description

This Anti-mouse IFN Alpha Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0369
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Mouse
Immunogen:	Recombinant Mouse IFN Alpha
Molecular Weight:	19.1 kDa (calculated) for mouse IFN- α
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Bovine IFN-A α : Weak; Caprine IFN-A: Moderate; Equine IFN α 1: Weak; Ovine IFN-A: Weak; Swine IFN α 1: Strong
Homology:	Mus musculus (house mouse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.