

Biotinylated Anti-equine IL-13 Polyclonal Antibody

KPAB0376

Description

This Biotinylated Anti-equine IL-13 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0376

Conjugated: Biotin

Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit

Reactivity: Equine

Immunogen: Recombinant Horse IL-13

Purification: Antigen-affinity Purification

Homology: Equus asinus (ass), Equus caballus (horse), Equus przewalskii (Przewalski's horse), Equus quagga (Plains Zebra)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.