

Anti-canine LIF Polyclonal Antibody

KPAB0377

Description

This Anti-canine LIF Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0377
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Canine
Immunogen:	Recombinant Canine LIF
Molecular Weight:	19.7 kDa (calculated) for canine LIF
Purification:	Antigen-affinity Purification
Formulation	1X PBS buffer containing 0.09% NaN ₃
Homology:	Canis lupus dingo (dingo), Canis lupus familiaris (dog), Vulpes lagopus (Arctic fox), Vulpes vulpes (red fox)
Cross-Reactivity:	Bovine LIF: Moderate; Caprine LIF: Moderate; Chicken LIF: Weak; Cynomolgus Monkey LIF: Moderate; Dolphin LIF: Moderate; Equine LIF: Weak; Feline LIF: Moderate; Guinea Pig LIF: Moderate; Human LIF: Moderate; Mouse LIF: Weak; Rabbit LIF: Moderate; Rat LIF: Moderate; Swine LIF: Moderate

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.