

## Biotinylated Anti-canine LIF Polyclonal Antibody

**KPAB0378**

### Description

---

This Biotinylated Anti-canine LIF Polyclonal Antibody has been validated for ELISA, ELISpot applications.

### Product Information

---

**SKU:** KPAB0378  
**Conjugated:** Biotin  
**Applications:** ELISA, ELISpot

### Antibody Data

---

**Host:** Rabbit  
**Reactivity:** Canine  
**Immunogen:** Recombinant Canine LIF  
**Purification:** Antigen-affinity Purification  
**Homology:** Canis lupus dingo (dingo), Canis lupus familiaris (dog), Vulpes lagopus (Arctic fox), Vulpes vulpes (red fox)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.