

Anti-rabbit LIF Polyclonal Antibody

KPAB0385

Description

This Anti-rabbit LIF Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0385
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Chicken
Reactivity:	Rabbit
Immunogen:	Recombinant Rabbit LIF
Molecular Weight:	19.7 kDa (calculated) for rabbit LIF
Purification:	Antigen-affinity Purification
Formulation	1X PBS buffer containing 0.09% NaN ₃
Homology:	Oryctolagus cuniculus (rabbit)
Cross-Reactivity:	Bovine LIF: Strong; Canine LIF: Strong; Caprine LIF: Strong; Chicken LIF: Weak; Cynomolgus Monkey LIF: Strong; Dolphin LIF: Strong; Equine LIF: Moderate; Feline LIF: Strong; Guinea Pig LIF: Strong; Human LIF: Strong; Mouse LIF: Strong; Rat LIF: Strong; Swine LIF: Strong

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.