

Anti-bovine EPO Polyclonal Antibody

KPAB0387

Description

This Anti-bovine EPO Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0387
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Bovine
Immunogen:	Recombinant Bovine Erythropoietin
Molecular Weight:	18.4 kDa (calculated) for bovine EPO
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Canine EPO: Strong; Equine EPO: Strong; Feline EPO: Strong; Mouse EPO: Strong; Rabbit EPO: Strong; Swine EPO: Strong
Homology:	Bison bison bison (American buffalo), Bos grunniens (domestic yak), Bos indicus (zebu), Bos mutus (wild yak), Bos taurus (cattle)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.