

Biotinylated Anti-bovine/caprine/ovine IL-5 Polyclonal Antibody

KPAB0394

Description

This Biotinylated Anti-bovine/caprine/ovine IL-5 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0394
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit
Reactivity: Bovine
Immunogen: Recombinant Bovine IL-5
Purification: Antigen-affinity Purification
Homology: Bison bison bison (American buffalo), Bos gaurus (gaur), Bos indicus (zebu), Bos taurus (cattle), Capra hircus (goat), Ovis aries (sheep)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.