

Anti-human CCL2 Polyclonal Antibody

KPAB0397

Description

This Anti-human CCL2 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0397
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Human
Immunogen:	Recombinant Human CCL2
Molecular Weight:	8.9 kDa (calculated) for human CCL2
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Bovine CCL2: Strong; Canine CCL2: Strong; Cynomolgus CCL2: Strong; Dolphin CCL2: Strong; Equine CCL2: Strong; Feline CCL2: Strong; Ferret CCL2: Moderate; Ovine CCL2: Moderate; Swine CCL2: Strong
Homology:	Gorilla gorilla gorilla (western lowland gorilla), Homo sapiens (human), Pan troglodytes (chimpanzee)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.