

## Biotinylated Anti-human CCL2 Polyclonal Antibody

**KPAB0398**

### Description

---

This Biotinylated Anti-human CCL2 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

### Product Information

---

**SKU:** KPAB0398

**Conjugated:** Biotin

**Applications:** ELISA, ELISpot

### Antibody Data

---

**Host:** Rabbit

**Reactivity:** Human

**Immunogen:** Recombinant Human CCL2

**Purification:** Antigen-affinity Purification

**Homology:** Gorilla gorilla gorilla (western lowland gorilla), Homo sapiens (human), Pan troglodytes (chimpanzee)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.