

## Anti-canine CXCL10 Polyclonal Antibody

**KPAB0401**

### Description

---

This Anti-canine CXCL10 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

### Product Information

---

<b>SKU:</b>	KPAB0401
<b>Concentration:</b>	1.0 mg/mL
<b>Conjugated:</b>	Unconjugated
<b>Applications:</b>	ELISA, ELISpot, FC, Neutralization, WB

### Antibody Data

---

<b>Host:</b>	Rabbit
<b>Reactivity:</b>	Canine
<b>Immunogen:</b>	Recombinant Canine CXCL10
<b>Molecular Weight:</b>	8.6 kDa (calculated) for canine CXCL10
<b>Purification:</b>	Antigen-affinity Purification
<b>Formulation:</b>	1X PBS buffer containing 0.09% NaN3
<b>Cross-Reactivity:</b>	Equine CXCL10: Weak; Feline CXCL10: Moderate; Mouse CXCL10: Moderate; Rabbit CXCL10: Moderate; Swine CXCL10: Strong
<b>Homology:</b>	Canis lupus dingo (dingo), Canis lupus familiaris (dog)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.