

Biotinylated Anti-canine CXCL10 Polyclonal Antibody

KPAB0402

Description

This Biotinylated Anti-canine CXCL10 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0402
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Rabbit
Reactivity:	Canine
Immunogen:	Recombinant Canine CXCL10
Molecular Weight:	8.6 kDa (calculated) for canine CXCL10
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN3
Cross-Reactivity:	Equine CXCL10: Weak; Feline CXCL10: Moderate; Mouse CXCL10: Moderate; Rabbit CXCL10: Moderate; Swine CXCL10: Strong
Homology:	Canis lupus dingo (dingo), Canis lupus familiaris (dog)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.