

Anti-mouse IFN gamma Polyclonal Antibody

KPAB0403

Description

This Anti-mouse IFN gamma Polyclonal Antibody has been validated for WB, ELISA applications.

Product Information

SKU: KPAB0403
Conjugated: Unconjugated
Applications: WB, ELISA

Antibody Data

Host: Rabbit
Reactivity: Mouse
Immunogen: Recombinant Mouse IFN gamma
Purification: Antigen-affinity Purification
Homology: Mus musculus (house mouse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.