

Anti-cynomolgus monkey IL-8 Polyclonal Antibody

KPAB0405

Description

This Anti-cynomolgus monkey IL-8 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

| | |
|-----------------------|--|
| SKU: | KPAB0405 |
| Concentration: | 1.0 mg/mL |
| Conjugated: | Unconjugated |
| Applications: | ELISA, ELISpot, FC, Neutralization, WB |

Antibody Data

| | |
|--------------------------|---|
| Host: | Rabbit |
| Reactivity: | Cynomolgus Monkey |
| Immunogen: | Recombinant Cynomolgus Monkey IL-8 |
| Molecular Weight: | 9.3 kDa (calculated) for cynomolgus monkey IL-8 |
| Purification: | Antigen-affinity Purification |
| Formulation: | 1X PBS buffer containing 0.09% NaN3 |
| Cross-Reactivity: | Human IL-8: Strong |
| Homology: | Macaca fascicularis (crab-eating macaque), Macaca mulatta (Rhesus monkey), Macaca nemestrina (pig-tailed macaque) |

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.