

Biotinylated Anti-cynomolgus monkey IL-8 Polyclonal Antibody

KPAB0406

Description

This Biotinylated Anti-cynomolgus monkey IL-8 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0406

Conjugated: Biotin

Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit

Reactivity: Cynomolgus Monkey

Immunogen: Recombinant Cynomolgus Monkey IL-8

Purification: Antigen-affinity Purification

Homology: Macaca fascicularis (crab-eating macaque), Macaca mulatta (Rhesus monkey), Macaca nemestrina (pig-tailed macaque)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.