

Biotinylated Anti-rabbit CXCL1 Polyclonal Antibody

KPAB0408

Description

This Biotinylated Anti-rabbit CXCL1 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0408
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Chicken
Reactivity: Rabbit
Immunogen: Recombinant Rabbit CXCL1
Purification: Antigen-affinity Purification
Homology: Oryctolagus cuniculus (rabbit)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.