

Biotinylated Anti-cynomolgus monkey CCL4 Polyclonal Antibody

KPAB0410

Description

This Biotinylated Anti-cynomolgus monkey CCL4 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0410
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Rabbit
Reactivity:	Cynomolgus Monkey
Immunogen:	Recombinant Cynomolgus Monkey CCL4
Purification:	Antigen-affinity Purification
Homology:	Macaca fascicularis (crab-eating macaque), Macaca mulatta (Rhesus monkey), Rhinopithecus roxellana (golden snub-nosed monkey)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.