

Biotinylated Anti-canine BAFF Polyclonal Antibody

KPAB0412

Description

This Biotinylated Anti-canine BAFF Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0412

Conjugated: Biotin

Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit

Reactivity: Canine

Immunogen: Recombinant Canine BAFF

Purification: Antigen-affinity Purification

Homology: Canis lupus dingo (dingo), Canis lupus familiaris (dog), Vulpes lagopus (Arctic fox), Vulpes vulpes (red fox)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.