

Biotinylated Anti-feline IFN alpha Polyclonal Antibody

KPAB0414

Description

This Biotinylated Anti-feline IFN alpha Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0414
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit
Reactivity: Feline
Immunogen: Recombinant Feline IFN alpha
Purification: Antigen-affinity Purification
Homology: Felis catus (domestic cat)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.