

Anti-human IL-1 beta Polyclonal Antibody

KPAB0417

Description

This Anti-human IL-1 beta Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0417
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Human
Immunogen:	Recombinant Human IL-1 beta
Molecular Weight:	17.4 kDa (calculated) for human IL-1 β
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Cynomolgus IL-1B: Weak; Rabbit IL-1B: Weak
Homology:	Homo sapiens (human), Pan paniscus (pygmy chimpanzee), Pan troglodytes (chimpanzee)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.