

Biotinylated Anti-human IL-1 beta Polyclonal Antibody

KPAB0418

Description

This Biotinylated Anti-human IL-1 beta Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0418

Conjugated: Biotin

Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit

Reactivity: Human

Immunogen: Recombinant Human IL-1 beta

Purification: Antigen-affinity Purification

Homology: Homo sapiens (human), Pan paniscus (pygmy chimpanzee), Pan troglodytes (chimpanzee)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.