

Biotinylated Anti-ovine IL-1 beta Polyclonal Antibody

KPAB0420

Description

This Biotinylated Anti-ovine IL-1 beta Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0420
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit
Reactivity: Ovine
Immunogen: Recombinant Sheep IL-1 beta
Purification: Antigen-affinity Purification
Homology: Ovis aries (sheep), Pseudois nayaur (bharal)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.