

Anti-ovine IL-4 Polyclonal Antibody

KPAB0421

Description

This Anti-ovine IL-4 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0421
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Ovine
Immunogen:	Recombinant Sheep IL-4
Molecular Weight:	12.6 kDa (calculated) for ovine IL-4
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Bovine IL-4: Weak; Caprine IL-4: Moderate
Homology:	Budorcas taxicolor (takin), Ovis aries (sheep)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.