

Anti-feline CXCL10 Polyclonal Antibody

KPAB0424

Description

This Anti-feline CXCL10 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0424
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Feline
Immunogen:	Recombinant Feline CXCL10
Molecular Weight:	8.7 kDa (calculated) for feline CXCL10
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Canine CXCL10: Weak; Swine CXCL10: Weak
Homology:	Acinonyx jubatus (cheetah), Felis catus (domestic cat), Lynx canadensis (Canada lynx), Lynx pardinus (Spanish lynx), Panthera leo (lion), Panthera pardus (leopard), Panthera tigris altaica (Amur tiger), Panthera uncia (snow leopard), Puma concolor (puma)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.