

Biotinylated Anti-dolphin CXCL10 Polyclonal Antibody

KPAB0425

Description

This Biotinylated Anti-dolphin CXCL10 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0425

Conjugated: Biotin

Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit

Reactivity: Dolphin

Immunogen: Recombinant Dolphin CXCL10

Purification: Antigen-affinity Purification

Homology: Delphinapterus leucas (beluga whale), Globicephala melas (long-finned pilot whale), Lagenorhynchus obliquidens (Pacific white-sided dolphin), Neophocaena asiaeorientalis asiaeorientalis (Yangtze finless porpoise), Orcinus orca (killer whale), Phocoena sinus (vaquita), Tursiops truncatus (common bottlenose dolphin)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.