

Biotinylated Anti-feline CXCL10 Polyclonal Antibody

KPAB0426

Description

This Biotinylated Anti-feline CXCL10 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0426

Conjugated: Biotin

Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit

Reactivity: Feline

Immunogen: Recombinant Feline CXCL10

Purification: Antigen-affinity Purification

Homology: Acinonyx jubatus (cheetah), Felis catus (domestic cat), Lynx canadensis (Canada lynx), Lynx pardinus (Spanish lynx), Panthera leo (lion), Panthera pardus (leopard), Panthera tigris altaica (Amur tiger), Panthera uncia (snow leopard), Puma concolor (puma)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.