

Anti-ferret CXCL10 Polyclonal Antibody

KPAB0427

Description

This Anti-ferret CXCL10 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0427
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Ferret
Immunogen:	Recombinant Ferret CXCL10
Molecular Weight:	8.5 kDa (calculated) for ferret/mink CXCL10
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Bovine CXCL10: Weak; Dolphin CXCL10: Moderate; Feline CXCL10: Moderate; Ovine CXCL10: Weak; Swine CXCL10: Moderate
Homology:	Mustela erminea (ermine), Mustela putorius furo (domestic ferret), Neogale vison (American mink)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.