

ACTB Monoclonal Antibody

MACO0090

Description

This ACTB Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	MACO0090
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Synonyms:	A26C1A antibody, A26C1B antibody, ACTB antibody, ACTB_HUMAN antibody, Actin beta antibody, Actin cytoplasmic 1 antibody, Actin, cytoplasmic 1, N-terminally processed antibody, Actx antibody, b actin antibody, b-actin antibody, Beta cytoskeletal actin antibody, Beta-actin antibody, BRWS1 antibody, E430023M04Rik antibody, Melanoma X actin antibody, MGC128179 antibody, PS1TP5 binding protein 1 antibody, PS1TP5BP1 antibody
Clone:	-
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Zebrafish

Antibody Data

Isotype:	IgG
Uniprot:	P60709
Host Species:	Mouse
Purification:	-
Immunogen:	Synthetic Peptide
Immunogen Species:	Homo sapiens (Human)
Buffer:	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:500-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western blot analysis of Zebrafish skeletal muscle, diluted at 1:3000.