

ACTB Monoclonal Antibody

MACO0180

Description

This ACTB Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	MACO0180
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Synonyms:	ACTB, BRWS1, PS1TP5BP1
Clone:	Monoclonal
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human,Rat,Mouse,Monkey,Dog

Antibody Data

Isotype:	-
Uniprot:	-
Host Species:	Mouse
Purification:	Antigen Affinity Purified
Immunogen:	-
Immunogen Species:	Homo sapiens (Human)
Buffer:	-
Form:	Mouse IgG1 in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:500-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-