

IgG2b Monoclonal Antibody -FITC Conjugated

MACO0323

Description

This IgG2b Monoclonal Antibody -FITC Conjugated is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	MACO0323
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Synonyms:	-
Clone:	Monoclonal
Applications:	ELISA IF FC
Conjugation:	FITC
Reactivity:	Human

Antibody Data

Isotype:	IgG2b
Uniprot:	-
Host Species:	Mouse
Purification:	-
Immunogen:	-
Immunogen Species:	Homo sapiens (Human)
Buffer:	-
Form:	0.2% BSA and 0.1% Proclin 300 Phosphate-buffered solution,PH7.2-7.4

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IF	1:100-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image

-

Description

-