

## IgG2b Monoclonal Antibody -PE Conjugated

MACO0324

### Description

---

This IgG2b Monoclonal Antibody -PE Conjugated is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	MACO0324
<b>Contents:</b>	100µl Bradford Reagent: 1 vial (2ml)
<b>Category:</b>	Monoclonal Antibody
<b>Synonyms:</b>	-
<b>Clone:</b>	Monoclonal
<b>Applications:</b>	ELISA IF FC
<b>Conjugation:</b>	PE
<b>Reactivity:</b>	-

### Antibody Data

---

<b>Isotype:</b>	IgG2b
<b>Uniprot:</b>	-
<b>Host Species:</b>	Mouse
<b>Purification:</b>	-
<b>Immunogen:</b>	-
<b>Immunogen Species:</b>	-
<b>Buffer:</b>	-
<b>Form:</b>	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA

## Preparation & Storage

---

**Storage:** Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IF	1:100-1:200

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data

---

Image	Description
-	-