

TNF Monoclonal Antibody

MACO0517

Description

This TNF Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	MACO0517
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Synonyms:	Cachectin,TNF-alpha,Tumor necrosis factor ligand superfamily member 2
Clone:	Monoclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Bovine

Antibody Data

Isotype:	IgG2b
Uniprot:	-
Host Species:	Mouse
Purification:	>95%,Protein G purified
Immunogen:	Recombinant Bovine TNF alpha protein
Immunogen Species:	Bos taurus (Bovine)
Buffer:	Preservative: 0.03% Proclin 300Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-