

Sdc1 Monoclonal Antibody

MACO0556

Description

This Sdc1 Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	MACO0556
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Synonyms:	Sdc1 antibody, Synd1 antibody, Syndecan-1 antibody, SYND1 antibody, CD antigen CD138 antibody
Clone:	Monoclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Rat

Antibody Data

Isotype:	IgG1
Uniprot:	P26260
Host Species:	Mouse
Purification:	>95%,protein G purified
Immunogen:	Recombinant Rat Syndecan-1 protein
Immunogen Species:	Rattus norvegicus (Rat)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-