

## Dehydroepiandrosterone(DHEA) Monoclonal Antibody

MACO0588

### Description

---

This Dehydroepiandrosterone(DHEA) Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	MACO0588
<b>Contents:</b>	50µg Bradford Reagent: 1 vial (2ml)
<b>Category:</b>	Small Molecules Antibodies
<b>Synonyms:</b>	DHEA, androstenolone, dehydroisoandrosterone, prasterone, diandrone, psicosterone, trans-dehydroandrosterone, diandron, 17-hormoforin
<b>Clone:</b>	Monoclonal
<b>Applications:</b>	For research use only
<b>Conjugation:</b>	Non-conjugated
<b>Reactivity:</b>	-

### Antibody Data

---

<b>Isotype:</b>	IgG1
<b>Uniprot:</b>	-
<b>Host Species:</b>	Mouse
<b>Purification:</b>	>95%, protein G purified
<b>Immunogen:</b>	Dehydroepiandrosterone-BSA conjugate
<b>Immunogen Species:</b>	-
<b>Buffer:</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Form:</b>	liquid

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

## Preparation & Storage

---

**Storage:** Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  
Store Bradford Reagent at Room Temperature for 1 Year.

<b>Recommended Dilutions:</b>	<b>Application</b>	<b>Recommended Dilution</b>
	-	-

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data

---

<b>Image</b>	<b>Description</b>
-	-