

Sulfadimidine Monoclonal Antibody

MACO0592

Description

This Sulfadimidine Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	MACO0592
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	Small Molecules Antibodies
Synonyms:	SM2, S-Mez, Sulfamethazine, 4,6-Dimethylsulfadiazine, 4-Amino-N-(4,6-dimethyl-2-pyrimidinyl)benzenesulfonamide
Clone:	Monoclonal
Applications:	For research use only
Conjugation:	Non-conjugated
Reactivity:	-

Antibody Data

Isotype:	IgG2a
Uniprot:	-
Host Species:	Mouse
Purification:	Caprylic Acid Ammonium Sulfate Precipitation purified
Immunogen:	Sulfadimidine-BSA conjugate
Immunogen Species:	-
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-