

PGC Antibody

PACO03211

Description

This PGC Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO03211
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	Polyclonal Antibody
Synonyms:	PEP1, Phosphoenolpyruvate carboxylase 1, PEPC 1, PEPCase 1, PPC
Clone:	-
Applications:	WB
Conjugation:	Non-conjugated
Reactivity:	Maize

Antibody Data

Isotype:	IgG
Uniprot:	P04711
Host Species:	Rabbit
Purification:	-
Immunogen:	Synthetic Peptide
Immunogen Species:	Maize
Buffer:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:2000-5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western blot analysis of Maize, diluted at 1:4000.