

NDUFC1 Antibody

PACO04407

Description

This NDUFC1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO04407
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	CI-KFYI antibody, Complex I KFYI subunit antibody, Complex I-KFYI antibody, KFYI antibody, mitochondrial antibody, NADH dehydrogenase (ubiquinone) 1, subcomplex unknown, 1, 6kDa antibody, NADH dehydrogenase [ubiquinone] 1 subunit C1 antibody, NADH dehydrogenase [ubiquinone] 1 subunit C1, mitochondrial antibody, NADH-ubiquinone oxidoreductase KFYI subunit antibody, NDUC1_HUMAN antibody, NDUFC1 antibody
Clone:	-
Applications:	IHC ELISA
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse

Antibody Data

Isotype:	IgG
Uniprot:	O43677
Host Species:	Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of Human NDUFC1.
Immunogen Species:	Homo sapiens (Human)
Buffer:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: Liquid

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:100-1:300
	ELISA	1:10000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-