

RALGAPA2 Antibody

PACO04798

Description

This RALGAPA2 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO04798
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	250 kDa substrate of Akt antibody, Akt substrate AS250 antibody, AS250 antibody, bA287B20.1 antibody, dJ1049G11 antibody, dJ1049G11.4 antibody, KIAA1272 antibody, p220 antibody, Ral GTPase-activating protein alpha subunit 2 antibody, Ral GTPase-activating protein subunit alpha-2 antibody, RALGAPA2 antibody, RapGAPalpha2 antibody, RGPA2_HUMAN antibody, RP11-470C13.2 antibody
Clone:	-
Applications:	WB IHC IF ELISA
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse

Antibody Data

Isotype:	IgG
Uniprot:	Q2PPJ7
Host Species:	Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of Human AS250.
Immunogen Species:	Homo sapiens (Human)
Buffer:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:500-1:2000
	IHC	1:100-1:300
	IF	1:200-1:1000
	ELISA	1:10000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-