

IKZF4 Antibody

PACO05030

Description

This IKZF4 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO05030
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	EOS antibody, Eos protein antibody, IKAROS family zinc finger 4 (Eos) antibody, IKAROS family zinc finger 4 antibody, Ikaros family zinc finger protein 4 antibody, IKZF4 antibody, IKZF4_HUMAN antibody, KIAA1782 antibody, Zinc finger protein Eos (IKAROS family zinc finger protein 4) antibody, Zinc finger protein Eos antibody, Zinc finger protein, subfamily 1A, 4 (Eos) antibody, Zinc finger protein, subfamily 1A, member 4 antibody, Zinc finger transcription factor Eos antibody, ZNFN1A4 antibody
Clone:	-
Applications:	WB ELISA
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse, Rat

Antibody Data

Isotype:	IgG
Uniprot:	Q9H2S9
Host Species:	Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the C-terminal region of Human Eos.
Immunogen Species:	Homo sapiens (Human)
Buffer:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: Liquid

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:500-1:2000
	ELISA	1:40000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-