

B3GALNT2 Antibody

PACO07960

Description

This B3GALNT2 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO07960
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti human PolyClonal Antibody
Synonyms:	3-GalNAc-T2 antibody, 3-N-acetylgalactosaminyltransferase 2 antibody, 3-N-acetylgalactosaminyltransferase II antibody, B3GalNAc T2 antibody, B3GalNAcT2 antibody, b3galnt2 antibody, B3GL2_HUMAN antibody, Beta 1,3 galactosaminyltransferase 2 antibody, Beta 1,3 N acetylgalactosaminyltransferase 2 antibody, Beta 1,3 N acetylgalactosaminyltransferase II antibody, Beta 1,3 N acetylgalactosaminyltransferase II MGC39558 antibody, Beta 3 GalNAc T2 antibody, Beta-1 antibody, beta-1,3-GalNAc-T2 antibody, Beta3GalNAcT2 antibody, MDDGA11 antibody, MGC39558 antibody, UDP GalNAc:betaGlcNAc beta 1,3 galactosa... Read more
Clone:	-
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse, Rat

Antibody Data

Isotype:	IgG
Uniprot:	Q8NCR0
Host Species:	Rabbit
Purification:	Antigen Affinity purified
Immunogen:	Human B3GALNT2
Immunogen Species:	Homo sapiens (Human)

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Buffer: PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Form: -

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-