

## AAMDC Antibody

PACO08073

### Description

---

This AAMDC Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                      |                                                                                                                                                       |
|----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>SKU:</b>          | PACO08073                                                                                                                                             |
| <b>Contents:</b>     | 100µl<br>Bradford Reagent: 1 vial (2ml)                                                                                                               |
| <b>Category:</b>     | Purified Rabbit Anti human PolyClonal Antibody                                                                                                        |
| <b>Synonyms:</b>     | AAMDC antibody, C11orf67 antibody, PTD015Mth938 domain-containing protein antibody, Adipogenesis associated Mth938 domain-containing protein antibody |
| <b>Clone:</b>        | -                                                                                                                                                     |
| <b>Applications:</b> | <b>ELISA</b> <b>IHC</b>                                                                                                                               |
| <b>Conjugation:</b>  | Non-conjugated                                                                                                                                        |
| <b>Reactivity:</b>   | Human, Mouse, Rat                                                                                                                                     |

### Antibody Data

---

|                           |                                                                                       |
|---------------------------|---------------------------------------------------------------------------------------|
| <b>Isotype:</b>           | IgG                                                                                   |
| <b>Uniprot:</b>           | Q9H7C9                                                                                |
| <b>Host Species:</b>      | Rabbit                                                                                |
| <b>Purification:</b>      | Antigen Affinity purified                                                             |
| <b>Immunogen:</b>         | Human C11orf67                                                                        |
| <b>Immunogen Species:</b> | Homo sapiens (Human)                                                                  |
| <b>Buffer:</b>            | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| <b>Form:</b>              | -                                                                                     |

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

## Preparation & Storage

---

**Storage:** Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  
Store Bradford Reagent at Room Temperature for 1 Year.

| Recommended Dilutions: | Application | Recommended Dilution |
|------------------------|-------------|----------------------|
|                        | -           | -                    |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data

---

| Image | Description |
|-------|-------------|
| -     | -           |