

PPDPF Antibody

PACO08101

Description

This PPDPF Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO08101
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti Human PolyClonal Antibody
Synonyms:	C20orf149 antibody, dJ697K14.9 antibody, Exdpf antibody, Exocrine differentiation and proliferation factor antibody, Pancreatic progenitor cell differentiation and proliferation factor antibody, Pancreatic progenitor cell differentiation and proliferation factor homolog antibody, PPDPF antibody, PPDPF_HUMAN antibody
Clone:	-
Applications:	ELISA WB IHC IF
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q9H3Y8
Host Species:	Rabbit
Purification:	Antigen Affinity purified
Immunogen:	Human C20orf149
Immunogen Species:	Homo sapiens (Human)
Buffer:	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Form:	-

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-