

DOC2B Antibody

PACO08907

Description

This DOC2B Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO08907
Contents:	DOC2B Antibody: 50 µl. Bradford Reagent: 1 vial (2ml)
Synonyms:	DOC2B antibody, Doc2-beta antibody, DOC2BL antibody, Double C2-like domain-containing protein beta antibody, Double C2-like domains beta antibody
Reactivity:	Human, Mouse, Rat
Applications:	ELISA, WB, IF
Form:	Liquid; PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.

Antibody Data

Immunogen:	Human DOC2B
Uniprot:	Q14184
Host Species:	Rabbit
Isotype:	IgG
Purification:	Antigen Affinity purified
Conjugated:	Unconjugated

Preparation & Storage

Storage: Store at -20°C or -80°C upon receipt. Avoid repeated freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit.
Visit www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.