

ENOPH1 Antibody

PACO09093

Description

This ENOPH1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO09093
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti human PolyClonal Antibody
Synonyms:	2 3 diketo 5 methylthio 1 phosphopentane phosphatase antibody, 2 antibody, 3-diketo-5-methylthio-1-phosphopentane phosphatase antibody, Acireductone synthase antibody, DKFZp586M0524 antibody, E 1 enzyme antibody, E1 antibody, Enolase phosphatase 1 antibody, Enolase phosphatase E1 antibody, Enolase-phosphatase E1 antibody, ENOPH_HUMAN antibody, Enoph1 antibody, FLJ12594 antibody, MASA homolog antibody, MST145 antibody, MSTP145 protein antibody
Clone:	-
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse, Rat

Antibody Data

Isotype:	IgG
Uniprot:	Q9UHY7
Host Species:	Rabbit
Purification:	Antigen Affinity purified
Immunogen:	Human ENOPH1
Immunogen Species:	Homo sapiens (Human)
Buffer:	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: -

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-