

GMEB2 Antibody

PACO09499

Description

This GMEB2 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO09499
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti human PolyClonal Antibody
Synonyms:	GMEB2 antibody, DNA binding protein p79PIF antibody, DNA-binding protein p79PIF antibody, Glucocorticoid modulatory element binding protein 2 antibody, Glucocorticoid modulatory element-binding protein 2 antibody, GMEB 2 antibody, GMEB-2 antibody, Gmeb2 antibody, GMEB2_HUMAN antibody, KIAA1269 antibody, P79PIF antibody, Parvovirus initiation factor p79 antibody, PIF p79 antibody, PIF79 antibody, RP4 663D7.1 antibody
Clone:	-
Applications:	ELISA WB IHC
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse, Rat

Antibody Data

Isotype:	IgG
Uniprot:	Q9UKD1
Host Species:	Rabbit
Purification:	Antigen Affinity Purified
Immunogen:	Human GMEB2
Immunogen Species:	Homo sapiens (Human)
Buffer:	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: -

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-