

IMP4 Antibody

PACO09948

Description

This IMP4 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO09948
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti human PolyClonal Antibody
Synonyms:	Brix domain containing protein 4 antibody, Brix domain-containing protein 4 antibody, BXDC4 antibody, D1Wsu40e antibody, IMP4 antibody, IMP4, U3 small nucleolar ribonucleoprotein, homolog (yeast) antibody, IMP4_HUMAN antibody, LRRGT00017 antibody, MGC19606 antibody, U3 small nucleolar ribonucleoprotein protein IMP4 antibody, U3 snoRNP protein 4 homolog antibody, U3 snoRNP protein IMP4 antibody
Clone:	-
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse, Rat

Antibody Data

Isotype:	IgG
Uniprot:	Q96G21
Host Species:	Rabbit
Purification:	Antigen Affinity purified
Immunogen:	Human IMP4
Immunogen Species:	Homo sapiens (Human)
Buffer:	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Form:	-

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-