

## INSC Antibody

PACO09967

### Description

---

This INSC Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                      |   |
|----------------------|---|
| <b>SKU:</b>          | PACO09967   |
| <b>Contents:</b>     | 100µl<br>Bradford Reagent: 1 vial (2ml)   |
| <b>Category:</b>     | Polyclonal Antibody   |
| <b>Synonyms:</b>     | Insc antibody, INSC_HUMAN antibody, Inscuteable homolog (Drosophila) antibody, Protein inscuteable homolog antibody |
| <b>Clone:</b>        | -   |
| <b>Applications:</b> | <a href="#">ELISA</a> <a href="#">WB</a> <a href="#">IHC</a> <a href="#">IF</a>                                     |
| <b>Conjugation:</b>  | Non-conjugated  |
| <b>Reactivity:</b>   | Human, Mouse, Rat   |

### Antibody Data

---

|                           |   |
|---------------------------|---|
| <b>Isotype:</b>           | IgG   |
| <b>Uniprot:</b>           | Q1MX18  |
| <b>Host Species:</b>      | Rabbit  |
| <b>Purification:</b>      | Antigen Affinity purified   |
| <b>Immunogen:</b>         | Human INSC  |
| <b>Immunogen Species:</b> | Homo sapiens (Human)  |
| <b>Buffer:</b>            | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| <b>Form:</b>              | Liquid  |

## Preparation & Storage

---

**Storage:** Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  
Store Bradford Reagent at Room Temperature for 1 Year.

| <b>Recommended Dilutions:</b> | <b>Application</b> | <b>Recommended Dilution</b> |
|-------------------------------|--------------------|-----------------------------|
|                               | -                  | -                           |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data

---

| <b>Image</b> | <b>Description</b> |
|--------------|--------------------|
| -            | -                  |