

LSM5 Antibody

PACO10303

Description

This LSM5 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO10303
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti human PolyClonal Antibody
Synonyms:	2010208O10Rik antibody, 2310034K10Rik antibody, FLJ12710 antibody, LSM5 antibody, LSM5 homolog antibody, LSM5 homolog, U6 small nuclear RNA associated (S. cerevisiae) antibody, LSM5_HUMAN antibody, U6 small nuclear RNA associated (S. cerevisiae) antibody, U6 snRNA-associated Sm-like protein LSm5 antibody, YER146W antibody
Clone:	-
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse, Rat

Antibody Data

Isotype:	IgG
Uniprot:	Q9Y4Y9
Host Species:	Rabbit
Purification:	Antigen Affinity Purified
Immunogen:	Human LSM5
Immunogen Species:	Homo sapiens (Human)
Buffer:	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Form:	-

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-